
HQ2D Modified Transfer Systems

Product Code	Description
HQ2D MOD MOT-N	Fully motorized transfer system meant for 2D transfers and also equipped with an integrated orientable needle that can pick up, displace and reorient a nanowire. The needle is connected to the software commanded mask holder.
HQ2D MOD MAN-N	Manual transfer system meant for 2D transfers and also equipped with an integrated orientable needle that can pick up, displace and reorient a nanowire. The needle is connected to the manually operated mask holder.
HQ2D MOD MAN-2XYZ	Manual transfer system with XY- and Z-stages for both the sample and the mask holder and an XY-stage for common displacement of sample and mask, so this system offers maximal freedom in the displacement of sample and/or mask. Can readily be equipped with stage motors.
HQ2D MOD MAN-2XYZ MROT	As HQ2D MOD MAN-2XYZ, with a motorized sample rotation stage.
HQ2D MOD MAN-BF	Standard manual transfer system with a pure bright field microscope without dark field illumination option.
